Nahita

Serie / Series / Série 601

DECLARACIÓN DE CONFORMIDAD CE
CE DECLARATION OF CONFORMITY
DÉCLARATION DE CONFORMITÉ CE

El fabricante | The manufacturer | Le fournisseur:
AUXILAB, S.L.

Declara que el equipo | Declare that the equipment | Declare que l'appareil:
BAÑOS TERMOSTÁTICOS | THERMOSTATIC WATER BATHS | BAINS THERMOSTATIQUES

Código | Code | Code:
50601003, 50601005, 50601012, 50601019

Serie | Series | Séries:
601

Cumple las siguientes directivas | Meet the following directives | Accomplit les directives suivantes:

2006/95/CE	Directiva sobre material eléctrico destinado a utilizarse con determinados límites de tensión (Directiva sobre Baja Tensión) Directive relating to electrical equipment designed for use within certain voltage limits (Low Voltage Directive) Directive relative aux matériels électriques à employer dans certaines limites de tension (Directive de basse tension)
2004/108/CE	Directiva de Compatibilidad Electromagnética (CEM) Directive for Electromagnetic Compatibility (EMC) Directive de Compatibilité Electromagnétique (CEM)

Cumple las siguientes Normas: | Meet the following Standards | Accomplit les normes suivantes:

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Thank you for choosing this equipment. We sincerely wish that you enjoy your Nahita digital thermostatic bath, series 601. We highly recommend looking after this equipment according to what is stated in this manual.

Nahita develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

Please bear in mind the following:

- ◆ This manual is inseparable from the Nahita thermostatic baths, so it should be available for all the users of these equipments.
- ◆ You should carefully handle the bath avoiding sudden movements, knocks, free fall of heavy / sharp objects on it.
- ◆ Never dismantle the different pieces of the bath to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the bath.
- ◆ To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting Nahita Technical Assistance Department by email to assistance@auxilab.es.
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the equipment's warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the warranty.
- ◆ Neither fuses nor accessories (including their loss), are covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Nahita Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Please do not forget filling the warranty certificate and send it before 15 days after the date of purchase.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.

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- ◆ Empty the tank using the emptying pump (except model 601/3).
- ◆ For cleaning the stainless steel tank use a cloth, never use scourers since they could scratch the tank. Take special care when cleaning the heating element and the temperature sensor.
- ◆ To avoid the appearance of calcareous deposits, it is recommended to clean the tank periodically with a mixture of water and vinegar (acetic acid) or with a solution of 10% of 28% chlorhydric acid. Rinse with clean water afterwards.

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- ◆ For cleaning the rest of the equipment, it is recommended to use soapy water without abrasive products.

6. STORAGE AND TRANSPORT

- ◆ Keep the equipment in a clean and dry environment with good ventilation and free from corrosive gases.
- ◆ During transport, prevent the equipment from wetting and crashing.

 **ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.**

 **INSTRUCTIONS ABOUT ENVIRONMENT PROTECTION**
Do not throw this equipment in the garbage when its life cycle had been completed; take it to a collection point for recycling electrical and electronic equipments. It does not contain dangerous or toxic elements for human health but it may cause damage to the environment. Materials may be recycled as indicated. If you recycle these materials, you are contributing to environment protection. Please contact the administration of your community to be advised on collection points.

 **ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.**

INDEX OF LANGUAGES

Spanish	2-11
English	12-20
French.....	21-29

INDEX OF CONTENTS

1. USES OF THE INSTRUMENT	13
2. DESCRIPTION	14
3. TECHNICAL SPECIFICATIONS	15
4. INSTALLATION / SETTING UP	15
5. MAINTENANCE AND CLEANING	19
6.STORAGE AND TRANSPORT	20
ANNEX I: CE CERTIFICATE	32

1. USES OF THE INSTRUMENT

Nahita digital thermostatic water baths have been designed to get a homogeneous temperature in samples and solutions. For a better performance, Nahita baths are thermally isolated between the tank and the external case, allowing a higher general performance, a better temperature homogenization and energy saving. The depth of the tank allows a bigger immersion and contact of samples with heating source, thus being appropriated to use in Biology, Chemistry, etc.


Thermostatic baths present an oven-painted external case that contains a stainless steel tank; inside this tank, at the lower part, there are the temperature sensor and the heating element, covered by a metallic rack to avoid direct contact between the recipients and the heating element.

The tank can be covered by the lid with holes and thanks to the reducing rings, get the appropriate diameter of the hole for any type of recipient. The diameter of the holes is approx. 11 cm and diameter of reducing rings is 8 and 5 cm.

Baths are equipped with a pump and drainpipe to be easily emptied (except model 601/3).


2. DESCRIPTION

- 1.1 Control panel
- 1.2 On/Off (Power): general on/off switch
- 1.3 On/Off (Outlet): emptying pump on/off switch
- 1.4 Drainpipe
- 1.5 Socket

Front panel

- 2.1 Set: parameter setting button
- 2.2 ▲: parameter increase button
- 2.3 ▼: parameter decrease button
- 2.4 ◀: digit selection button
- 2.5 PV display (Parameter Value)
- 2.6 SV display (Set value)
- 2.7 RUN indicator
- 2.8 OUT indicator: heating element working
- 2.9 AT indicator: autoadjustment
- 2.10 ST indicator: set time
- 2.11 ALM indicator: alarm

- ◆ Be careful when touching the equipment. The metallic tank can reach temperatures up to 100°C. Once the equipment has been turned off it can remain hot for several minutes
- ◆ Pay attention to the set temperature; never leave unattended the equipment during operation and with heating function on.
- ◆ Temperature must always be set to at least 50°C lower than the ignition point of the media used.
- ◆ Pay attention to the safety and resistance to high temperatures of the container used. Process pathogenic materials only in closed vessels.
- ◆ During operation wear the necessary personal protections to avoid the risk of possible damages as burns caused by splashing and evaporation of liquids or intoxication caused by release of toxic or combustible gases.
- ◆ Before operation, always check the instrument and accessories are not damaged. Do not use damaged components. Safe operation is only guaranteed with the accessories described in the "Listing of items" chapter. Accessories must be securely attached to the device and cannot come off by themselves. Always disconnect the plug before fitting accessories.

5. MAINTENANCE AND CLEANING

- ◆ Proper maintenance helps to keep instruments working in a good state and lengthen its lifetime.

Note: all the processes of use mention above will not have any value unless it is kept a continued and careful maintenance.

- ◆ Please follow the processes of use of this manual.
This manual should be available for all users of this equipment.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ The equipment is supplied with a standard Schuko wire; it has to be plugged to an earth connection and the socket should be handy and ready to unplug the equipment in case of emergency.
- ◆ If the equipment is not going to be used during a long period of time; make sure it is unplugged.
- ◆ Never try to repair the equipment by yourself, since warranty will be lost and damages in the equipment, users (burns, hurts,...), electrical installation or nearby equipments can be caused.
- ◆ In the event of breakdown please contact your distributor to overhaul through Nahita Technical Assistance Department

Cleaning

Note: before cleaning the bath, disconnect it from the main supply and let the heating element cool down.

Alarm

In case real temperature exceeds in more than 10°C the set temperature, the alarm will sound; at the same time real temperature will be displayed on PV display (2.5), command dHAL will be displayed on SV display (2.6) and ALM indicator (2.11) will flick (red).

◆ Turn the equipment off by pressing On/Off switch (Power) (1.2) to stop the alarm sound.

Emptying

All baths except model 601/3 are equipped with a pump and drainpipe (1.4) for an easy emptying.

Note: before emptying the bath make sure the On/Off switch (Power) (1.2) is in off position (O); otherwise the heating element could be working when the tank is empty.

◆ Press On/Off switch (Outlet) (1.3) and the pump will start emptying the tank through the drainpipe (1.4).

Auto-adjust

The auto-adjust function allows the adjustment of the internal parameters of the regulator to a certain working temperature. For this proceed as follows:

- ◆ Select the working temperature as explained on “*Temperature setting*” of this manual.
- ◆ Keep pressed ▼ button (2.3) until the AT indicator (2.9) starts flicking; auto-adjust has begun.
- ◆ Auto-adjust will finish when AT indicator (2.9) lights off.
- ◆ After the auto-adjust the equipment will operate perfectly.

Security

- ◆ Read the user’s manual before operating with the equipment. This equipment must be used only by qualified staff.
- ◆ Connect the device to an earthed power supply to ensure safety of machine and experiment; the voltage indicated on the equipment must correspond to that of the power supply used.
- ◆ The instrument can only be disconnected from the main power supply by pulling out the connector plug, not the cable.
- ◆ This equipment is forbidden to use with inflammable, explosive, poisonous and strong corrosive solutions.
- ◆ Make sure the equipment is placed and installed in a plain, horizontal and stable surface.
- ◆ Respect the minimum safety distances between devices, between the device and the wall and above the assembly (min.200 mm).
- ◆ Keep the equipment away from high magnetic field.
- ◆ The bath must be only used with the tank full of water; if the heating element operates without being covered by water it will be damaged.

3. TECHNICAL SPECIFICATIONS

Nahita digital thermostatic water baths have been designed to get a homogeneous temperature in samples and solutions. For a better performance, Nahita baths are thermally isolated between the tank and the external case, allowing a higher general performance, a better temperature homogenization and energy saving. The depth of the tank allows a bigger immersion and contact of samples with heating source, thus being appropriated to use in Biology, Chemistry, etc.

Thermostatic baths present an oven-painted external case that contains a stainless steel tank; inside this tank, at the lower part, there are the temperature sensor and the heating element, covered by a metallic rack to avoid direct contact between the recipients and the heating element.

The tank can be covered by the lid with holes and thanks to the reducing rings, get the appropriate diameter of the hole for any type of recipient. The diameter of the holes is approx. 11 cm and diameter of reducing rings is 8 and 5 cm.

Baths are equipped with a pump and drainpipe to be easily emptied (except model 601/3).

<i>Model</i>	601/3	601/5	601/12	601/19
<i>Reference</i>	50601003	50601005	50601012	50601019
<i>Capacity</i>	2.5 L	5 L	12 L	19 L
<i>Useful dimensions (cm)</i>	15x13x10	15x30x10	32x30x10	50x30x10
<i>External dimensions (cm)</i>	18x17x21	33x18x21	36x33x22	54x34x22
<i>Holes</i>	1	2	4	6
<i>Max. temperature</i>	100 °C			
<i>Resolution</i>	1 °C			
<i>Accuracy</i>	± 1 °C			
<i>Consumption</i>	300 W	500 W	1000 W	1000 W
<i>Fuse</i>	5 A	5 A	8 A	10 A
<i>Power</i>	220 V ± 10%; 50/60 ± 1 Hz			

4. INSTALLATION / SETTING UP**Preliminary inspection**

- ◆ Unwrap the bath and take off the polyspan protection in which it comes fitted. Take off the involving plastic and without connecting the bath to the power supply, make sure that the equipment does not present any damage because of the shipment. In case it presents any damage, tell it immediately to your transport agent or dealer so that they can make the claims in the correct time limit.
- ◆ Please keep the original wrapping; you will always need it for returns enclosed with all the accessories supplied.

- ◆ Please check that all the accessories are enclosed with the equipment:
 - Bath
 - Power cable
 - Lid with holes
 - Reducing rings with lid (depending on model)
 - Rack to cover heating element
 - User's manual

Note: we will only accept the equipment return provided it comes in its original wrapping.

Installation

- ◆ Before using the instrument, you should become familiar with its components and basic basis, as well as functions of its controls.

READ CAREFULLY THIS MANUAL BEFORE OPERATING WITH THIS EQUIPMENT, TO GET THE MAXIMUM AND A LONGER LIFE.

- ◆ Place the bath on a horizontal, plane and stable surface appropriate for the weight of the equipment.
- ◆ Make sure the switches of the equipment are in Off (0) position.
- ◆ Make sure the required operating voltage and power supply voltage match; ensure the socket is earthed reliably.
- ◆ Put the metallic rack so as it covers the heating element and fill the tank up to 2/3 its capacity.

Note: It is recommended to use decalcified water. DO NOT use deionized or demineralized water since they have a high rusting power. To prolong bath lifetime and avoid rusting, dissolve 1 g of CO₃Na₂ (sodium carbonate) per liter of water used in the bath.

Setting up

BEFORE OPERATING WITH THE BATH, CHECK THE LEVEL OF WATER IN THE TANK; IF IT IS NECESSARY ADD MORE WATER. BEAR IN MIND THAT IF THE HEATING ELEMENT IS WORKING WITHOUT BEING COVERED BY WATER IT WILL BE DAMAGED.

- ◆ Press the On/Off switch (1.2).
- ◆ PV display (2.5) will show real temperature of the bath, while SV display (2.6) will show the set temperature value.
- ◆ After its use do not forget turning off the bath by pressing the On/Off switch (1.2) and unplug it from the power supply.

- ◆ If the equipment is not going to be used during a long period of time, empty the tank with the pump and drainpipe (1.4).

Temperature setting

- ◆ The maximum working temperature is 100° C; proceed as follows to set working temperature:
 - ◆ Press Set (2.1); command "Sv" will be displayed on PV display (2.5) and temperature set value on SV display (2.6).
 - ◆ Use ▲ (2.2) and ▼ (2.3) buttons to select the desired temperature value (show non SV display). Use ◀ (2.4) button to move from one digit to the next.
 - ◆ Press Set (2.1). The heating element will start working to reach the set temperature; PV display will show the real temperature of the water, while SV display will show the set temperature.
 - ◆ While the heating element is working, RUN indicator (2.7) will be on (green) and OUT indicator (2.8) will flick (Green); once the set temperature is reached then OUT indicator (2.8) will light off.

Note: to get an homogeneous temperature in all the volume of water, stir regularly the water in the tank.

Note: to assure a constant and homogeneous temperature and prevent water evaporation from the bath, place the lid and the appropriate reducing rings to the size of the recipient to be warmed.

Time setting

By default, when turning the equipment on the time will be off, so as the equipment will be operating until it is turned off by pressing On/Off switch (1.2). In case it is necessary to work during periods of time, proceed as follows:

- ◆ Press ▲ button (2.1). Command "St" will be displayed on PV display (2.5) and by default "Off" will be displayed on SV display (2.6).
- ◆ Use ▲ (2.2), ▼ (2.3) and ◀ (2.4) buttons to set the working time value; time value can be set up to 9999 min.
- ◆ Press Set (2.1) to save the time value; the heating element will star working.
- ◆ During operation RUN (2.7) and ST (2.10) indicators will be on (green) and OUT indicator (2.8) will flick (green).
- ◆ Time count will start as soon as set temperature is reached.

Note: if working time is long, keep an eye on water level to avoid the heating element working without being covered by water.

When working time is finished the equipment and heating element will stop heating (all indicators will be off). No alarm will sound when time finishes.

- ◆ To restart the bath, keep pressed ◀ button (2.4) until RUN indicator (2.7) lights on. A new working temperature and time can be set (time will be again Off by default).

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958 412 886

629 636 705

<http://www.cromakit.es/>

Calle Tucumán 8 Nave B 18200 Maracena (Granada)